


CALL

sound alerting device

Tongeber

unité d'alerte sonore

dispositivo de alerta sonora



Ⓞ user guide

Ⓞ Bedienungsanleitung

Ⓞ instructions

Ⓞ instrucciones

))))) **PhonicEar Logia™**

other resources online at www.phonicearlogia.com

CALL function

CALL is the PUZZLE system's sound alerting device. A sound signal is sent by CALL's speaker when an alarm from one of the PUZZLE system's alerting transmitters (DETECT and/or OBSERVER) is received. CALL can receive up to seven different alarms, apart from the alarm for battery change. Each alarm has its own unique sound, which makes it easier to distinguish between the different alarm signals. The individual sounds are selected on the alerting transmitter (DETECT and/or OBSERVER).

CALL features

- ① 12V AC/DC connection
To connect 12volt power transformer

- ② volume control
To adjust volume level

- ③ channel tuner
To select personal code


- ④ tone control
To select frequency

- ⑤ indicator light
Shows that power is on


- ⑥ speaker
Sounds signal received from other PUZZLE units

- ⑦ power transformer
Supplies power to CALL unit



- ⑧ screwdriver
Used to set personal code and frequency



back




front




CALL accessories (available separately)


9 table stand
Allows CALL unit to be placed on a table or other flat surface



10 wall bracket
Allows CALL unit to be mounted on the wall



11 back plate
To be used with the table stand and the wall bracket



CALL instructions

1a positioning CALL unit



We recommend that you install CALL in the room you use most often.

CALL can either be hung on a wall or used with its table stand (sold separately) and remain portable.


1b positioning CALL unit



CALL can be affixed to the wall with Velcro tape, a screw, or a wall bracket (sold separately).

CALL instructions

2 connecting transformer



1 The cable from the enclosed transformer is inserted into the socket marked **12V AC/DC** on the back of CALL.

2 The transformer should then be connected to the wall socket, and the wall contact switched on.

When CALL is installed correctly, the red light on the front of CALL should come on.

3 select personal code



The personal code is selected using the tuner marked **Channel** on the back of CALL.

NOTE: *The PUZZLE system's different transmitters and receivers should be set to the same personal code.*

NOTE: *Use the screwdriver included with CALL to turn dial to desired setting.*

CALL instructions

4 select the frequency

Tone Tuner	Frequency (Hz)
0	125
1	187
2	250
3	300
4	350
5	400
6	500
7	625
8	750
9	875
A	1000
B	1400
C	1800
D	2000
E	3000
F	4000

The sound signal should be set to the frequency that the user hears best.

If you already know that frequency, follow the scheme on the left.

If you do not know the frequency, try listening to each of the 16 frequencies until you find the one you hear best.

The frequency is selected by using the tuner marked **Tone** on the back of CALL.

NOTE: Use the screwdriver included with CALL to turn dial to desired setting.

5 set the volume



CALL's volume should be adjusted to best suit the user's needs.

The volume control is marked **Volume** on the back of CALL.

Turn the control clockwise to increase the volume, and counter-clockwise to decrease it.

troubleshooting

Should your CALL unit not function satisfactorily, please check the list below:

The red light on the front of CALL unit does not light

- Check to see that the wall switch it turned on.
- Switch off the wall contact for 10 seconds and then switch it back on again.

The red light on the front of CALL lights, but you cannot hear a signal when one of the PUZZLE system's alerting transmitters is activated

- Check that the volume is set correctly. See [5].
- Check that the personal code on CALL and the PUZZLE systems' alerting transmitters are the same. See [3].

CALL emits an alarm without apparent reason

- Change the personal code using the tuner marked **Channel** on the back of CALL. See [3].

If you still have trouble after trying the ideas above please contact your local authorized PhonicEar Logia dealer/installer (see information on back cover).

test and approval

CALL is approved and CE marked in accordance with the R&TTE directive.

warranty

A one-year limited warranty is issued to the original consumer of this product(s) by PhonicEar Logia.

time period of warranty

This warranty will go into effect upon the date of original purchase of this equipment, and will remain in effect for one calendar year from that date as long as the instrument remains the property of the original consumer.

what is covered by this warranty

Any electronic component which, because of workmanship, manufacturing, or design defects, fails to function properly under normal use during the life of this warranty will be replaced or repaired at no charge for parts or labor when returned to the factory service center. Transportation is paid by the customer. If it is determined that repair is not feasible, the entire unit may be replaced with an equivalent unit upon mutual agreement of the manufacturer and customer.

what is not covered by this warranty

This one-year limited warranty does not apply to:

- 1** Malfunctions resulting from abuse, neglect, or accident.
- 2** Peripheral accessories as itemized within the product brochure, when such items are returned after 90 days from the original purchase.
- 3** Alkaline batteries (if applicable).
- 4** Instruments connected, installed, used, or adjusted in a manner contrary to instructions provided by the manufacturer.
- 5** Consequential damages and damages resulting from delay or loss of this instrument. The exclusive remedy under this warranty is strictly limited to repair or replacement as herein provided.
- 6** Products damaged in transit unless investigated by the shipper and returned to the warrantor with the investigation report.


PhonicEar Logia reserves the right to make changes in the design or construction of any of its instruments at any time without incurring any obligation to make any changes whatsoever on units previously purchased. This warranty is in lieu of all other expressed warranties. All expressed and implied warranties will terminate upon the expiration of this written warranty. No representative or person is authorized to represent or assume for us any liability in connection with the sale or use of our products other than as set forth above.

The above warranty does not affect any legal rights that you might have under applicable national legislation governing sale of consumer goods. Your distributor may have issued a warranty that goes beyond the clauses of this limited warranty. Please consult your distributor for further information.

what to do if you need service

If you require service under the terms of this warranty, carefully package the instrument to prevent damage and return it to your distributor (see information on back cover). Please provide your receipt, a detailed description of the problem, your full name, billing address, and telephone number.

Your local dealer:



europe: +45 3917 7101 • **germany (kostenlose Rufnummer):** 0800/1844736

u.s.: 800.227.0735, then press 7 • **canada:** 800.263.8700

other regions: Call your local PhonicEar Logia dealer • **website:** www.phonicearlogia.com