


WAKE

vibrating alarm and alerting device

Weck- und Alarmsystem

alarme et réveil matin vibrants

alarma de vibración y dispositivo de alerta



(GB) user guide

(D) Bedienungsanleitung

(F) instructions

(E) instrucciones

)))))) **PhonicEar Logia™**

other resources online at www.phonicearlogia.com

WAKE function

WAKE is the PUZZLE system's vibrating alarm and alerting device. By using a powerful external vibrator WAKE makes the user aware of the alarm clock ringing.

WAKE also alerts you when an alarm from one of the PUZZLE system's alerting transmitters is received (DETECT or OBSERVER alerting transmitter). WAKE can receive up to seven different alarms apart from the alarm signaling a low battery. Each alarm has its own unique vibration pattern, which makes it easier to distinguish between the different alarms. The individual vibration patterns are selected on the alerting transmitter (DETECT and/or OBSERVER).

WAKE features


① 12V AC connection
To connect 12volt power transformer

② vibrator connection
To connect vibrator

③ channel tuner
To set personal code

④ indicator light
Shows that power is on

⑤ microphone sensitivity tuner
To adjust microphone sensitivity



back




front




bottom

WAKE accessories (included with WAKE system)

- ⑥ vibrator
Vibrates to alert user that alarm clock is ringing or when an alarm from one of the PUZZLE system's alerting transmitters is received



- ⑦ power transformer
Supplies power to WAKE unit




- ⑧ screwdriver
Used to set personal code and adjust microphone sensitivity tuner



WAKE accessories (not included, available separately)

- ⑨ battery backup
To supply power to WAKE unit if power outage should occur



- ⑩ WAKE alarm clock
Photo not available


WAKE instructions

1 position WAKE unit



WAKE is placed on the bedside table, with its microphone as close as possible to the alarm clock.

2 connect transformer



1 Plug transformer cord into **12V AC/DC** connection on back of WAKE.

2 Plug transformer into wall socket and turn on the contact. Red light on front of WAKE should go on.

WAKE instructions

3 connect vibrator



Plug vibrator cable into the **Vib** connection on back of WAKE.

4 position vibrator



Place vibrator under pillow.

WAKE instructions

5 select personal code



Use the **Channel** tuner on the back of WAKE to set personal code.

NOTE: *The PUZZLE system's different alerting transmitters and receivers should be set to the same personal code.*

NOTE: *Use the screwdriver included with WAKE to turn dial to desired settings.*

6a adjust microphone sensitivity



Change sensitivity by adjusting the control that is visible through the oblong holes on the bottom of WAKE.

Turn the control counterclockwise to reduce sensitivity and clockwise to increase sensitivity.


The control should only be moved very slightly between tests.

NOTE: *The WAKE's sensitivity is set at the factory and should not be changed unless there are false alarms or problems with the alarm clock not activating WAKE.*

NOTE: *Use the screwdriver included with WAKE to turn dial to desired settings.*


WAKE instructions

6b) adjust microphone sensitivity



Sensitivity is correctly adjusted once the red light on the front of WAKE blinks or when the vibrator vibrates after the alarm clock rings.

7) installing optional battery backup




The battery backup ensures that WAKE will work even if there is a power loss.

- 1 Plug battery backup cord into **12V AC/DC** connection on back of WAKE
- 2 Plug transformer cord into battery backup connection
- 3 Plug transformer into wall socket and turn on the contact

Red light on front of WAKE unit and green light on end of battery backup unit should both go on.

8 battery backup function



1 If a power outage occurs, the green light at one end of the battery backup unit will go out.

2 The red light on the front of WAKE will continue to shine.

When fully charged the battery backup unit can supply WAKE with enough power to run for 24 hours on standby and for approximately three minutes when the vibrator is active.

The battery backup unit recharges automatically when the power supply is re-established so it is important that the wall contact be on all the time.

troubleshooting

If your WAKE vibrating alarm and alerting device does not function satisfactorily, try the solutions listed below:

The red light on the front of WAKE does not go on.

- Check that the wall contact is switched on and working.
- Turn the wall contact off for 10 seconds, then turn it back on and look for the light.

The red light on the front of WAKE goes on, but the vibrator does not vibrate when a PUZZLE system alerting transmitter is activated.

- Check that the personal code on WAKE and the PUZZLE system's alerting transmitters are the same. See [5].

The red light on the front of WAKE lights and the vibrator vibrates when one of PUZZLE system's transmitters is activated, but not when the alarm clock rings.

- Check that the microphone sensitivity is correctly adjusted. See [6a] and [6b].

Vibrator vibrates without apparent reason.

- Check that the microphone's sensitivity is correctly adjusted. See [6a] and [6b].
- Change the personal code by using the tuner marked **Channel** on the back of WAKE. See [5].

If you still have trouble after trying the ideas above please contact your local authorized PhonicEar Logia dealer/installer.

test and approval

WAKE is approved and CE marked in accordance with the R&TTE directive.

warranty

A one-year limited warranty is issued to the original consumer of this product(s) by PhonicEar Logia.

time period of warranty

This warranty will go into effect upon the date of original purchase of this equipment, and will remain in effect for one calendar year from that date as long as the instrument remains the property of the original consumer.

what is covered by this warranty

Any electronic component which, because of workmanship, manufacturing, or design defects, fails to function properly under normal use during the life of this warranty will be replaced or repaired at no charge for parts or labor when returned to the factory service center. Transportation is paid by the customer. If it is determined that repair is not feasible, the entire unit may be replaced with an equivalent unit upon mutual agreement of the manufacturer and customer.

what is not covered by this warranty

This one-year limited warranty does not apply to:

- 1** Malfunctions resulting from abuse, neglect, or accident.
- 2** Peripheral accessories as itemized within the product brochure, when such items are returned after 90 days from the original purchase.
- 3** Alkaline batteries (if applicable).
- 4** Instruments connected, installed, used, or adjusted in a manner contrary to instructions provided by the manufacturer.
- 5** Consequential damages and damages resulting from delay or loss of this instrument. The exclusive remedy under this warranty is strictly limited to repair or replacement as herein provided.
- 6** Products damaged in transit unless investigated by the shipper and returned to the warrantor with the investigation report.


PhonicEar Logia reserves the right to make changes in the design or construction of any of its instruments at any time without incurring any obligation to make any changes whatsoever on units previously purchased. This warranty is in lieu of all other expressed warranties. All expressed and implied warranties will terminate upon the expiration of this written warranty. No representative or person is authorized to represent or assume for us any liability in connection with the sale or use of our products other than as set forth above.

The above warranty does not affect any legal rights that you might have under applicable national legislation governing sale of consumer goods. Your distributor may have issued a warranty that goes beyond the clauses of this limited warranty. Please consult your distributor for further information.

what to do if you need service

If you require service under the terms of this warranty, carefully package the instrument to prevent damage and return it to your distributor (see information on back cover). Please provide your receipt, a detailed description of the problem, your full name, billing address, and telephone number.

Your local dealer:



europa: +45 3917 7101 • **germany (kostenlose Rufnummer):** 0800/1844736

u.s.: 800.227.0735, then press 7 • **canada:** 800.263.8700

other regions: Call your local PhonicEar Logia dealer • **website:** www.phonicearlogia.com