



Phonic Ear is a Danish owned company in William Demant Holding A/S. The group has more than 100 years experience in hearing and sound and comprises a number of international companies.



A simple concept for all alarms and call signals in your home

When you lose some of your hearing, things which people with normal hearing take for granted can suddenly cause problems.

Puzzle makes everyday life easier for anyone with impaired hearing. Our assistive listening devices ensure that you can:

- hear anyone ringing your doorbell
- hear your phone ringing
- hear your alarm clock
- hear your baby crying
- hear your smoke alarm

Discreet, flexible design

Puzzle product design is contemporary and nice, and our products are made in sizes to fit in perfectly with any home. They are flexible and can be placed and moved, however you like.

Your local distributor:



© 2014 Phonic Ear A/S

Produced by:



Kongebakken 9 ▪ 2765 Smørum ▪ Denmark
 tel: +45 3917 7101 ▪ fax: +45 3927 7900
 mail@phonicear.com ▪ www.phonicear.com

921-37-410-00
 Rev. E 0314

Hear all alarms and call signals in your home



Hear all alarms and call signals in your home

Alerting transmitters

Choose one or more alarm transmitters depending on the alarm/call you want to hear.

DETECT



Wireless alerting transmitter that provides notification of various alarms/call signals, or when your smoke alarm is going off.

1

OBSERVER



Wireless baby and/or call alarm that notifies you when your baby is crying or a person in bed needs help.

2

DETECT DOOR



Wireless doorbell that notifies you when somebody is at the door.

3

Alerting receivers

Choose one or more alarm receiver depending on how you want to be notified of alarms/call signals: with sound, light or vibration.

FLASH



Visual light alerting device that uses a powerful light signal to notify you of alarms/call signals.

4

CALL



Sound receiver that uses a powerful and frequency-adapted audible signal to notify you of alarms/call signals.

5

POCKETVIB



Portable wireless vibration receiver that uses powerful vibrations to notify you of alarms/call signals.

6

WAKE



Wake-up system that uses a powerful vibrator to notify you of alarms/call signals when you are sleeping.

7

